

## **Renafem+ EMS +RF**

- ✓ 3 in 1 Technologies
- ✓ +35% Muscle Building
- ✓ 15 Tesla High Intensity

# **RenaSculpt**

Model: FE60



WhatsApp: +8618533632296



Email: [winkonlaser@winkonlaser.com](mailto:winkonlaser@winkonlaser.com)

Web: [www.winkonlaser.com](http://www.winkonlaser.com) / [www.beautemax.com](http://www.beautemax.com)

Add: Yanjiao Yinhe Bay Youth Development Base, Beijing, China




## Workout Smarter Not Harder

Experience muscle gain in just 30 minutes. Our machine stimulates muscle growth by promoting higher frequency contractions, surpassing traditional exercise. Four treatment sessions equal 12-16 weeks of training, with microcurrent and micro RF technology enhancing blood circulation.


### EMS + Micro-RF

### Renafem



Large Applicator for flat areas for example, abdomen

Small Applicator for curve areas like arms



### High-Intensity Focused Electro-Magnetic Technology

Utilizing microcurrent and micro radio frequency technology, our machine simultaneously emits high-frequency magnetic contraction signals and releases microcurrents to directly destroy fat cell tissues. The damaged fat cells are naturally eliminated, reducing fat layer thickness, and sculpting a beautiful silhouette.





### Pelvic Chair Cushion (Optional)

- The active regeneration of collagen in pelvic floor tissues, strengthens pelvic floor muscles, and lifts buttocks
- Improve strength, neuromuscular control, and sexual intimacy

### • Comprehensive Renting System

This system is specifically designed to manage and control the machine's usage time, number of uses, and operation dates. Primarily intended for our distributors, it facilitates the process of renting out machines to their clients. By enabling precise tracking and control of machine usage, the Comprehensive Renting System ensures optimal and fair rental practices, contributing to efficient business operations.

### • First-Class ABS Matte Shell



The machine's exterior is encased in a first-class ABS (Acrylonitrile Butadiene Styrene) shell with a matte finish. This high-grade ABS material is known for its robustness and durability, ensuring that our machine is built to last while also offering a sleek, professional appearance.



## FAQ

### Q: How long does it take to see results?

A: The timeframe can vary as it depends on individual factors. However, the majority of users typically start observing results within 3 to 5 weeks of consistent use.





### Q: What areas can be treated with Renasculpt?

- Neck / Shoulders
- Triceps / Biceps / Forearms
- Glutes / Abs
- Quadriceps / Hamstrings / Calves


### Q: Is this technology safe?

A: Yes, this technology has been approved by the FDA and is non-invasive, making it a safe option.



### Q: What might a patient experience during the treatment?

During the treatment, patients might experience intense muscle contractions and vibrations, accompanied by a mild sensation of electrical stimulation.







## Specification

Model	FE60
Technology	Renafem+ EMS+ Micro-current RF
Handle Piece	Large flat applicators *2 & Small curve aplicators*2 Pelvic Floor Pad *1 (Optional)
Magnetic Vibration Intensity	1-15 Tesla
Renafem Frequency	1-200HZ
Treatment Mode	Combined Mode& Alternate Mode
Input Voltage	AC 110V-220V
Output Power	5000W
Pulse Duration	350μs
Fuse	20A
Cooling System	Air Cooling
Size of Flight Shipping Case	69*57*112cm
Gross /Net Weight	90/71KG
Working Mode	Automatic Mode and Manual Mode